



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human ADAM-like Protein Decysin-1/ADAMDEC1 Protein(C-6His)

产品货号: TP05093

第三版

别名:ADAM DEC1;a disintegrin and metalloproteinase domain-like protein decysin-1;ADAM-like protein decysin-1

描述:ADAM DEC1 protein is expressed highly in the small intestine and appendix, moderately in lymph node, mucosal lining of the colon, thymus, spleen and very weakly in the bone marrow. ADAM DEC1 is induced during DC maturation and up-regulated in response to T-cell signals. It may play an important role in the control of the immune response and during pregnancy.

配方:PBS

分子量:30.3 kDa

序列:Lys206-Glu470

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.