

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Adaptin Ear-Binding Coat-Associated Protein 2/NECAP2 Protein(N, C-6His)

产品货号: TP05094

第三版

别名:Adaptin ear-binding coat-associated protein 2; NECAP endocytosis associated protein 2; NECAP2

描述:NECAP2 belongs to the NECAP family. The WXXF motifs mediate binding of accessory proteins to the ear-domain of AP-1, GGAs and AP-2 through hydrophobic interactions. Adaptin ear-binding coat-associated protein 2 can interact with AP1G1 and AP2A1 components of the adapter protein complex AP-1 and AP-2. It also interacts with the GAE domain proteins GGA1, GGA2 and GGA3.

配方:PBS

分子量:31.5 kDa

序列:Met 1-Phe263

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.