



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Polymeric Immunoglobulin Receptor/PIgR Protein(C-Fc)

产品货号: TP05102

第三版

别名:Polymeric Immunoglobulin Receptor; PIgR; Poly-Ig Receptor; Hepatocellular Carcinoma-Associated Protein TB6; PIGR

描述:PIGR binds polymeric IgA and IgM at the basolateral surface of epithelial cells. The complex is then transported across the cell to be secreted at the apical surface. During this process a cleavage occurs that separates the extracellular (known as the secretory component) from the transmembrane segment.

配方:PBS

分子量:95.0 kDa

序列:Lys19-Arg638

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.