

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Limbic system-associated membrane protein Protein

产品货号: TP05103

第三版

别名: Limbic system-associated membrane protein; LSAMP; IgLON family member 3

描述: Limbic system-associated membrane protein is also known as LSAMP, IgLON family member 3. In humans, it is encoded by the LSAMP gene. It belongs to the immunoglobulin superfamily and contains 3 Ig-like C2-type domains. Limbic system-associated membrane protein mediates selective neuronal growth and axon targeting. It contributes to the guidance of developing axons and remodeling of mature circuits in the limbic system. It is also essential for normal growth of the hippocampal mossy fiber projection

配方: PBS

分子量: 32.8 kDa

序列: Val29-Asn315

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.