



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Roundabout homolog 4 Protein

产品货号: TP05107

第三版

别名: Roundabout homolog 4; Magic roundabout; ROBO4; UNQ421/PRO3674

描述: Roundabout homolog 4 (ROBO4) contains 2 fibronectin type-III domains, 2 Ig-like C2-type (immunoglobulin-like) domains and belongs to the immunoglobulin superfamily. The protein is a receptor for slit proteins, at least for SLIT2, and seems to be involved in angiogenesis and vascular patterning. ROBO4 may mediate the inhibition of primary endothelial cell migration by Slit proteins.

配方: PBS

分子量: 74.6 kDa

序列: Gln28-Arg467

纯度: > 95% by HPLC

浓度:

内毒素: < 1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.