

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Hyaluronic Acid-binding Protein

产品货号: TP05110

第三版

别名: Complement Component 1 Q Subcomponent-Binding Protein Mitochondrial; ASF/SF2-Associated Protein p32; Glycoprotein gC1qBP; C1qBP; Hyaluronan-Binding Protein 1; Mitochondrial Matrix Protein p32; gC1q-R Protein; p33; C1QBP; GC1QBP; HABP1; SF2P32

描述: Complement Component 1Q Subcomponent-Binding Protein (C1QBP) is a nucleus protein that belongs to the MAM33 family. C1QBP is known to bind to the globular heads of C1q molecules and inhibit C1 activation. Mitochondrial C1QBP is a critical mediator of p14ARF-induced apoptosis. C1QBP functions as a chemotactic factor for immature dendritic cells, and migration is mediated through ligation of both C1QBP and cC1qR/CR. C1QBP overexpression successfully blocks mRNA accumulation from the adenovirus major late transcription unit (MLTU) and stimulates RNA polymerase II carboxy-terminal domain phosphorylation in virus-infected cells.

配方: PBS

分子量: 24.9 kDa

序列: Leu74-Gln282

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.