

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse S100 Calcium Binding Protein A15A/S100A15A Protein

产品货号: TP05122

第三版

别名:S100 calcium-binding protein A15A; Protein S100-A15A; Protein S100-A7A; S100 calcium-binding protein A7A; S100a15a

描述:Members of the S100 protein family are involved in calcium- or zinc-dependent cellular functions and regulate immune-mechanisms, cell proliferation and differentiation. Some S100 members have been established as tumor markers because they are dysregulated during carcinogenesis. Psoriasin (S100A7) and koebnerisin (S100A15) are highly homologous proteins that have been first described in psoriasis, which is characterized by disturbed epidermal maturation and chronic inflammation. Several studies showed that the coexpression of the hS100A7 and hS100A15 in psoriasis suggests that both proteins participate in keratinocyte maturation, proliferation and/or skin inflammation.

配方:PBS

分子量:12.9 kDa

序列:Met1-Tyr108

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.