
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse SLAM Family Member 9/SLAMF9/CD2F-10 Protein(C-6His)

产品货号: TP05137

第三版

别名:SLAM family member 9;Slamf9

描述:Mouse SLAM family member 9(Slamf9) is a single-pass type I membrane protein. The SLAMF9 gene encodes a member of the signaling lymphocytic activation molecule family. The encoded protein is a cell surface molecule that consists of two extracellular immunoglobulin domains, a transmembrane domain and a short cytoplasmic tail that lacks the signal transduction motifs found in other family members. The slamf9 protein may play a role in the immune response.

配方:PBS

分子量:24.7 kDa

序列:Phe18-Leu230

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.