

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human LAMP1/CD107a Protein(C-Fc)

产品货号: TP05146

第三版

**别名:**Lysosome-Associated Membrane Glycoprotein 1; LAMP-1; Lysosome-Associated Membrane Protein 1; CD107 Antigen-Like Family Member A; CD107a; LAMP1

**描述:**Lysosome-Associated Membrane Glycoprotein 1 (LAMP1) is a single-pass type I membrane protein belonging to the LAMP family. LAMP1 is expressed largely in the endosome-lysosome membranes of cells. It shuttles between lysosomes, endosomes, and the plasma membrane. LAMP1 functions to present carbohydrate ligands to selectins and it has also been implicated in tumor cell metastasis. It has been proposed LAMP1 can be used as a therapeutic agent for certain cancers, as well as a marker for lysosomal storage disorders and degranulation on lymphocytes such as CD8+ and NK cells. Cell surface LAMP1 and LAMP2 have been shown to promote adhesion of human peripheral blood mononuclear cells(PBMC) to vascular endothelium, therefore they are possibly involved in the adhesion of PBMCs to the site of inflammation.

**配方:**PBS

**分子量:**65.5 kDa

**序列:**Ala29-Met382

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.