



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human E-Selectin//ELAM1/CD62E Protein(C-6His)

产品货号: TP05151

第三版

**别名:**E-selectin (SELE); CD62 antigen-like family member E; Endothelial leukocyte adhesion molecule 1; Leukocyte-endothelial cell adhesion molecule 2; CD62E and ELAM1

**描述:**SELE is expressed on the surface of endothelial cells and mediates the interaction of leukocytes and platelets with endothelial cells during an inflammatory response. SELE functions as a cell-surface glycoprotein and has a role in immunoadhesion. In addition, SELE may also have a role in capillary morphogenesis.

**配方:**PBS

**分子量:**58.6 kDa

**序列:**Trp22-Pro556

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.