

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Eukaryotic Translation Initiation Factor 5A-2/EIF5A2 Protein(N-6His)

产品货号: TP05154

第三版

别名:Eukaryotic translation initiation factor 5A-2; Eukaryotic initiation factor 5A isoform 2; EIF5A2

描述:EIF5A2 is a member of the eIF-5A family. It is a 153 amino acids protein that in humans is encoded by the EIF5A2 gene. EIF5A2 is expressed in ovarian and colorectal cancer cell lines and highly expressed in testis. It has an important function at the level of mRNA turnover and mediates effects of polyamines on neuronal process extension and survival. It plays an important role in brain development and function, and in skeletal muscle stem cell differentiation.

配方:PBS

分子量:18.9 kDa

序列:Met 1-Lys153

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.