

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human BOC Protein/BOC Protein(C-6His)

产品货号: TP05162

第三版

别名:BOC protein;Boc homolog (Mouse); isoform CRA_b;Brother of CDO;

描述:Brother of CDO is a protein that in humans is encoded by the BOC gene. CDON and BOC are cell surface receptors of the immunoglobulin (Ig)/fibronectin type III repeat family involved in myogenic differentiation. CDON and BOC are coexpressed during development, form complexes with each other in a cis fashion, and are related to each other in their ectodomains, but each has a unique long cytoplasmic tail.

配方:PBS

分子量:14.2 kDa

序列:Asp31-Ser157

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.