



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Leydig Insulin-Like 3/INSL3 Protein(C-6His)

产品货号: TP05167

第三版

别名:Insulin-like 3; Leydig insulin-like peptide; Relaxin-like factor; INSL3; RLF; RLNL

描述:Insulin-like 3 is a protein that in humans is encoded by the INSL3 gene. It is a secreted protein that belongs to the insulin family. It is expressed in prenatal and postnatal Leydig cells and found as well in the corpus luteum, trophoblast, fetal membranes and breast. It may act as a hormone to regulate growth and differentiation of gubernaculum, and thus mediating intra-abdominal testicular descent. It is a ligand for LGR8 receptor.

配方:PBS

分子量:13.4 kDa

序列:Leu21-Tyr131

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.