

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Regenerating Islet-Derived Protein 4/RELP Protein(C-6His)

产品货号: TP05174

第三版

别名:Regenerating islet-derived protein 4; Gastrointestinal secretory protein; REG-like protein; Regenerating islet-derived protein IV; GISP; RELP; REG4

描述:REG4 is a secreted contains one C-type lectin domain, high expressed in the gastrointestinal tract, including jejunum, ileum, appendix, pancreas and small intestine. REG4 can be up-regulated by mucosal injury from active Crohn ' s disease or ulcerative colitis. In the acid environment, REG4 can maintain carbohydrate recognition activity. REG4 may be involved in inflammatory and metaplastic response of the gastrointestinal epithelium.

配方:PBS

分子量:19.2 kDa

序列:Asp23-Pro158

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.