

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Retinol-Binding Protein 3 Protein(N-6His)

产品货号: TP05181

第三版

别名:Retinol-binding protein 3; Interphotoreceptor retinoid-binding protein; IRBP; Interstitial retinol-binding protein; RBP3

**描述:**Retinol-binding proteins (RBP) are a family of proteins with diverse functions. They are carrier proteins that bind retinol. Retinol and retinoic acid play crucial roles in the modulation of gene expression and overall development of an embryo. However, deficit or excess of either one of these substances can cause early embryo mortality or developmental malformations. Regulation of transport and metabolism of retinol necessary for a successful pregnancy is accomplished via RBP. Retinol binding proteins have been identified within the uterus, embryo, and extraembryonic tissue of the bovine, ovine, and porcine, clearly indicating that RBP plays a role in proper retinol exposure to the embryo and successful transport at the maternal-fetal interface.

**配方:**PBS

**分子量:**35.2 kDa

**序列:**Thr321-Leu630

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.