

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Secretogranin-2/SCG2 Protein(C-6His)

产品货号: TP05186

第三版

别名:Secretogranin-2;Chromogranin-C;Secretogranin II; SgII.

描述:Secretogranin-2 is also known as Chromogranin-C,Secretogranin II, SgII. In humans, it is encoded by the SCG2 gene.It belongs to the chromogranin/secretogranin protein family. Secretogranin-2 is a neuroendocrine secretory granule protein, which is the precursor for biologically active peptides. It derived secretoneurin was distributed with strong immunoreactivity in the somata of pelvic ganglion neurons,72% of which also contained tyrosine hydroxylase,as well as in nerve terminals in the muscular layer and the lamina propria of the vas deferens.

配方:PBS

分子量:68.8 kDa

序列:Ser29-Met617

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.