

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human C-C Motif Chemokine 11/CCL11/Eotaxin Protein

产品货号: TP05202

第三版

别名:Eotaxin; C-C Motif Chemokine 11; Eosinophil Chemotactic Protein; Small-Inducible Cytokine A11; CCL11; SCYA11

描述:C-C Motif Chemokine 11 (CCL11) is a secreted protein that belongs to the intercrine beta (chemokine CC) family. In response to the presence of allergens, CCL11 selectively recruits eosinophils, a prominent feature of allergic inflammatory reactions. The effects of CCL11 are mediated by its binding to a G-protein-linked receptor known as a chemokine receptor. Chemokine receptors for CCL11 include CCR2, CCR3 and CCR5. However, it has been found that CCL11 has high degree selectivity for its receptor, such that they are inactive on neutrophils and monocytes, which do not express CCR3.

配方:PBS

分子量:8.6 kDa

序列:Gly24-Pro97

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.