



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Galectin-Related Protein/LGALSL Protein(C-6His)

产品货号: TP05203

第三版

别名:Galectin-Related Protein; Lectin Galactoside-Binding-Like Protein; LGALSL; GRP; HSPC159

描述:Galectin-Related Protein (LGALSL) is a 172 amino acid protein that contains one Galectin domain. LGALSL does not appear to bind carbohydrates or lactose as the critical residues required for binding are not conserved. LGALSL may play a significant role in stimulating smooth muscle growth in developing alveolar wall vessels and the development of pulmonary capillaries.

配方:PBS

分子量:20.0 kDa

序列:Ala2-Ser172

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.