

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse Fc receptor-like protein 1 Protein

产品货号: TP05208

第三版

别名:Fc receptor-like protein 1; FcR-like protein 1; FcRL1; BXMAS1-like protein 1; mBXMH1; Fc receptor homolog 1; FcRH1 ; moFcRH1 ; IFGP family protein 1; mIFGP1; CD307a

描述:Mouse Fc receptor-like protein 1 (FCRL1) is a single-pass type I membrane protein. It is expressed in B-cells at the various stages of differentiation. Mouse FCRL1 consists of a 203 amino acid (aa) extracellular domain (ECD) with two Ig-like domains, a 21 aa transmembrane segment, and a 103 aa cytoplasmic domain with two immunotyrosine activation motifs (ITAMs). FCRL1 is likely not a receptor for immunoglobulin. It may function as an activating coreceptor in B-cells and B-cells differentiation.

配方:PBS

分子量:22.5 kDa

序列:Ala17-Ser219

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.