

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Mouse C-X-C Motif Chemokine 16/CXCL16/SR-PSOX Protein(C-6His)

产品货号: TP05248

第三版

**别名:**C-X-C motif chemokine 16;Scavenger receptor for phosphatidylserine;oxidized low density lipoprotein; Small-inducible cytokine B16; Transmembrane chemokine CXCL16; SR-PSOX; Zmynd15

**描述:**CXCL16 is a single-pass type I membrane protein, which consists of 246 amino acids, CXCL16 induces a strong chemotatic response and calcium mobilization. CXCL16 acts as a scavenger receptor on macrophages, which specially binds to oxidized low density lipoprotein. CXCL16 may involves in pathophysiology such as atherogenesis. Soluble CXCL16 may play an important role in liver metastases through the induction of epithelial-mesenchymal transition.

**配方:**PBS

**分子量:**20.1 kDa

**序列:**Asn27-Trp201

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.