

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Nucleolar Protein 3/NOL3 Protein(N-GST)

产品货号: TP05278

第三版

别名:Nucleolar Protein 3; Apoptosis Repressor With CARD; Muscle-Enriched Cytoplasmic Protein; Myp; Nucleolar Protein of 30 kDa; Nop30; NOL3; ARC; NOP

描述:Nucleolar Protein 3 is encoded by NOL3 gene; multiple transcript variants encoding different isoforms have been found for this gene. So far, Nucleolar protein 3 has show to have two Isoforms. Isoform 1 may be involved in RNA splicing.Isoform 2 may inhibit apoptosis.It has been shown to down-regulate the enzyme activities of caspase 2, caspase 8 and tumor protein p53.

配方:PBS

分子量:48.9 kDa

序列:Met 1-Ser208

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.