

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Ornithine Decarboxylase/ODC1 Protein(C-6His, N-T7 tag)

产品货号: TP05287

第三版

别名:Ornithine decarboxylase; ODC1;

描述:Ornithine decarboxylase (ODC1) is an enzyme which belongs to the Orn/Lys/Arg decarboxylase class-II family. ODC1 is the rate-limiting enzyme of the polyamine biosynthesis pathway which catalyzes ornithine to putrescine. The activity level for the enzyme varies in response to growth-promoting stimuli and exhibits a high turnover rate in comparison to other mammalian proteins.

配方:PBS

分子量:53.5 kDa

序列:Met 1-Val461

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.