



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Corneodesmosin/CDSN Protein(C-6His)

产品货号: TP05293

第三版

别名:Corneodesmosin; S Protein; CDSN

描述:Corneodesmosin is a protein that is encoded by the CDSN gene in humans. This gene encodes a protein found in corneodesmosomes, which localize to human epidermis and other cornified squamous epithelia. During maturation of the cornified layers, the protein undergoes a series of cleavages, which are thought to be required for desquamation.

配方:PBS

分子量:49.3 kDa

序列:Lys33-Pro528

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.