



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Pancreatic Lipase-Related Protein 2/PLRP2 Protein(C-6His)

产品货号: TP05300

第三版

别名:Pancreatic Lipase-Related Protein 2; PL-RP2; Galactolipase; PNLIPRP2; PLRP2

描述:Pancreatic Lipase-Related Protein 2 (PNLIPRP2) is a secreted protein that belongs to the Lipase family of AB hydrolase superfamily. PNLIPRP2 is a lipase with broad substrate specificity that contains one PLAT domain. PNLIPRP2 is expressed mainly in pancreas. PNLIPRP2 can hydrolyze both phospholipids and galactolipids. PNLIPRP2 acts preferentially on monoglycerides, phospholipids and galactolipids. PNLIPRP2 also contributes to milk fat hydrolysis.

配方:PBS

分子量:51.2 kDa

序列:Lys18-Cys469

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.