

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Tumor Necrosis Factor Receptor II/TNFRSF1B/CD120b Protein(288-461, C-6His)

产品货号: TP05332

第三版

别名:Tumor necrosis factor receptor superfamily member 1B;Tumor necrosis factor receptor 2;TNF-R2;Tumor necrosis factor receptor type II;TNF-RII;TNFR-II;p75;p80 TNF-alpha receptor

描述:Tumor necrosis factor receptor superfamily member 1B(TNFRSF1B) is expressed by the gene TNFRSF1B. The soluble form is produced from the membrane form by proteolytic processing. It can bind to TRAF2, and interacts with BMX. It can act as the receptor with high affinity for TNFSF2/TNF-alpha and approximately 5-fold lower affinity for homotrimeric TNFSF1/lymphotoxin-alpha. The TRAF1/TRAF2 complex recruits the apoptotic suppressors BIRC2 and BIRC3 to TNFRSF1B/TNFR2. This receptor mediates most of the metabolic effects of TNF-alpha.

配方:PBS

分子量:19.0 kDa

序列:Lys288-Ser461

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.