

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human 40S Ribosomal Protein S19/RPS19 Protein

产品货号: TP05334

第三版

别名:40S Ribosomal Protein S19; RPS19

描述:40S Ribosomal Protein S19 (RPS19) is a ribosomal protein that belongs to the ribosomal protein S19e family. RPS19 is located in the nucleoli, and higher level expression is seen in colon carcinoma tissue than normal colon tissue. It required for pre-rRNA processing and maturation of 40S ribosomal subunits. RPS19 plays a role in many biological processes, such as endocrine pancreas development, erythrocyte differentiation, mRNA metabolic process. Defects in RPS19 are the cause of Diamond-Blackfan anemia type 1 (DBA1), which is a form of Diamond-Blackfan anemia, a congenital non-regenerative hypoplastic anemia that usually presents early in infancy. Diamond-Blackfan anemia is characterized by a moderate to severe macrocytic anemia, erythroblastopenia, and an increased risk of malignancy.

配方:PBS

分子量:16.1 kDa

序列:Pro2-His145

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.