

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Golgi Membrane Protein 1 Protein

产品货号: TP05339

第三版

别名: Golgi Membrane Protein 1; Golgi Membrane Protein GP73; Golgi Phosphoprotein 2; GOLM1; C9orf155; GOLPH2

描述: Golgi Membrane Protein 1 (GOLM1) belongs to the GOLM1/CASC4 family, GOLM1 is a single-pass type II membrane protein and can be found in many tissues. GOLM1 is overexpressed in prostate cancer and lung adenocarcinoma tissue. GOLM1 can be up-regulated in response to viral infection. GOLM1 is induced by the E1A adenoviral protein. GOLM1 plays a key role in the sorting and modification of proteins exported from the endoplasmic reticulum.

配方: PBS

分子量: 42.6 kDa

序列: Ser35-Leu400

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.