

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Interferon Lambda-2/IL-28A Protein(C-6His)

产品货号: TP05341

第三版

别名: Interferon lambda-2; IFN-lambda-2; Cytokine Zcyto20; Interleukin-28A; IL-28A; IFNL2; IL28A; ZCYTO20

描述: IL-28A (Interferon- $\lambda$ 2/IFN- $\lambda$ 2), IL-28B/IFN- $\lambda$ 3, and IL-29/IFN- $\lambda$ 1 are type III interferons which are distantly related to IL-10 family and type I IFN family cytokines. Mature human IL-28A is an approximately 22-25 kDa protein that shares 66% amino acid (aa) sequence identity with mouse and rat IL-28A and shows cross-species activity. It shares 96% and 70% aa sequence identity with human IL-28B and IL-29, respectively. IL-28A promotes the Th1 polarization of dendritic cells in the airway and inhibits Th2 and Th17 mediated inflammation. IL-28A additionally exhibits anti-tumor activity, in part by enhancing IL-12 dependent anti-tumor CTL responses in vivo. In contrast, it is up-regulated in invasive bladder cancer where it promotes tumor cell migration.

配方: PBS

分子量: 20.6 kDa

序列: Val26-Val200

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.