

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Ubiquitin-Conjugating Enzyme E2 B/UBE2B/HR6B Protein(C-6His)

产品货号: TP05355

第三版

别名: Ubiquitin-Conjugating Enzyme E2 B; RAD6 Homolog B; HR6B; hHR6B; Ubiquitin Carrier Protein B; Ubiquitin-Conjugating Enzyme E2-17 kDa; Ubiquitin-Protein Ligase B; UBE2B; RAD6B

描述: Ubiquitin-Conjugating Enzyme E2 B (UBE2B) is a member of the ubiquitin-conjugating enzyme family. UBE2B accepts ubiquitin from the E1 complex and catalyzes its covalent attachment to other proteins. It is shown that UBE2B interacts with RAD18, UBR2, and WAC. UBE2B is required for post-replicative DNA damage repair. Additional, UBE2B plays a role in sepsis-induced muscle protein proteolysis and cancer-induced cachexia.

配方: PBS

分子量: 18.3 kDa

序列: Met 1-Ser152

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.