



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Cytoglobin/CYGB Protein(C-6His)

产品货号: TP05358

第三版

别名:Cytoglobin; Histoglobin; HGb; Stellate Cell Activation-Associated Protein; CYGB; STAP

描述:Cytoglobin is a ubiquitously globin protein that belongs to the globin family. The highest expressed in heart, stomach, bladder and small intestine. CYGB acts a protector under conditions of oxidative stress. CYGB may be involved in intracellular oxygen storage or transfer, modulates oxygen and nitric oxide metabolism or scavenging free radicals within a cell.

配方:PBS

分子量:22.4 kDa

序列:Met 1-Pro190

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.