



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Deoxyhypusine Synthase/DHS Protein(C-6His)

产品货号: TP05373

第三版

别名:Deoxyhypusine Synthase; DHS; DHPS; DS

描述:Human Deoxyhypusine Synthase (DHS) is vital for the first step of hypusine biosynthesis. DHS catalyzes the NAD-dependent oxidative cleavage of spermidine, the subsequent transfer of the butylamine moiety of spermidine to the epsilon-amino group of a specific lysine residue of the eIF-5A precursor protein to form the intermediate deoxyhypusine residue.

配方:PBS

分子量:42.0 kDa

序列:Met 1-Asp369

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.