

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Dickkopf-Related Protein 3/DKK3 Protein(C-6His)

产品货号: TP05384

第三版

别名: Dickkopf-Related Protein 3; Dickkopf-3; Dkk-3; hDkk-3; DKK3; REIC

描述: Dickkopf-related protein 3 (DKK3) belongs to the DKK protein family including Dkk-1, 2, 3 and -4. DKK3 is a 350 amino acid secreted glycoprotein which is comprised of an N-terminal signal peptide and 2 conserved cysteine-rich domains that are separated by a 12 amino acid linker region. Dkk-3 also have one prokineticin domain. DKK3 is involved in embryonic development through its inhibition of the WNT signaling pathway. The Dkk family also includes Soggy, which is homologous to Dkk-3 but not to the other family members. Soggy has not been shown to inhibit Wnt signaling, and its role in the pathway is unclear.

配方: PBS

分子量: 37.2 kDa

序列: Ala22-Ile350

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.