



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Polyubiquitin-B/UBB Protein(N-GST, 6His)

产品货号: TP05388

第三版

别名:Polyubiquitin-B;UBB

描述:Polyubiquitin-B(UBB) is one of Ubiquitins. Ubiquitin is one of the most conserved proteins known in eukaryotic organisms. Ubiquitin is required for ATP-dependent, non-lysosomal intracellular protein degradation of abnormal proteins and normal proteins with a rapid turnover. Ubiquitin is covalently bound to proteins to be degraded, and presumably labels these proteins for degradation.

配方:PBS

分子量:9.4 kDa

序列:Met 153-Gly228

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.