



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human VIP36-Like Protein/LMAN2L Protein(C-6His)

产品货号: TP05401

第三版

别名:VIP36-like protein;Lectin mannose-binding 2-like;LMAN2-like protein;VIPL

描述:VIP36-like protein (LMAN2L) is a single-pass type I membrane protein and contains 1 L-type lectin-like domain. It is highly expressed in skeletal muscle and kidney, and its intermediate expression levels in heart, liver and placenta, low levels in brain, thymus, spleen, small intestine and lung. LMAN2L may be involved in the regulation of export from the endoplasmic reticulum of a subset of glycoproteins. It also may function as a regulator of ERGIC-53.

配方:PBS

分子量:34.4 kDa

序列:Ser19-Ala313

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.