

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Prion-Like Protein Doppel/PRND Protein(C-6His)

产品货号: TP05404

第三版

别名: Prion-like protein doppel; PrPLP; Prion protein 2; PRND; DPL

描述: Prion-like protein doppel is a 152 amino acids protein that belongs to the prion family. The protein encoded by this gene is found on chromosome 20, approximately 20 kbp downstream of the gene encoding cellular prion protein, to which it is biochemically and structurally similar. It is a membrane glycosylphosphatidylinositol-anchored glycoprotein that is found predominantly in testis. Mutations in this gene may lead to neurological disorders.

配方: PBS

分子量: 15.5 kDa

序列: Arg27-Gly152

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.