



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human V-Set and Ig Domain-Containing Protein 8/VSIG8 Protein(C-Fc)

产品货号: TP05411

第三版

别名:V-set and immunoglobulin domain-containing protein 8;VSIG8;C1orf204

描述:V-set and immunoglobulin domain-containing protein 8(VSIG8) is a single-pass type I membrane protein. The human VSIG8 cDNA encodes 414 amino acids (aa) including a 21 aa signal sequence, a 242 aa extracellular domain (ECD) containing 2 Ig-like V-type (immunoglobulin-like) domains, a 21 aa transmembrane domain and a 130 aa cytoplasmic domain. The function of VSIG8 is not clear.

配方:PBS

分子量:54.2 kDa

序列:Val22-Gly263

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.