

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human V-Set and Transmembrane Domain-Containing 1/VSTM1 Protein(C-6His)

产品货号: TP05413

第三版

别名:V-set and transmembrane domain-containing protein 1;SIRL-1;Signal inhibitory receptor on leukocytes-1;VSTM1

描述:V-set and transmembrane domain-containing protein 1 is a single-pass membrane protein. VSTM1 Contains 1 Ig-like V-type domain, in humans is encoded by the VSTM1 gene. It expressed on myeloid (neutrophils, eosinophils and monocytes) but not on lymphoid cells. It behaves as a cytokine, promoting IL17A secretion by CD4+ T-cells, and differentiation and activation of IL17 producing helper T-cells (TH17). Inhibitory immune receptor involved in the regulation of phagocytes.

配方:PBS

分子量:14.6 kDa

序列:Tyr17-Thr135

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.