
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Annexin A7/ANXA7 Protein

产品货号: TP05414

第三版

别名:Annexin A7; Annexin VII; Annexin-7; Synexin; ANXA7; ANX7; SNX

描述:Annexin A7 (ANXA7) is a member of the annexin family of calcium-dependent phospholipid binding proteins. Annexin A7 has a unique, highly hydrophobic N-terminal domain and a conserved C-terminal region. The C-terminal region is composed of alternating hydrophobic and hydrophilic segments. Annexin A7 is a calcium/phospholipid-binding protein with diverse properties including voltage-sensitive calcium channel activity and promotes membrane fusion and is also involved in exocytosis.

配方:PBS

分子量:50.3 kDa

序列:Met 1-Gln466

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.