

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Annexin A8/ANXA8 Protein

产品货号: TP05417

第三版

别名:Annexin A8; Annexin VIII; Annexin-8; Vascular Anticoagulant-Beta; VAC-Beta; ANXA8; ANX8

描述:Annexin A8 (ANXA8) belongs to a family of calcium-dependent membrane and phospholipid binding proteins. Annexin A8 contains 4 annexin repeats separated by linking sequences of variable lengths. A pair of annexin repeats may form one binding site for calcium and phospholipid. ANXA8 is preferentially expressed in acute promyelocytic leukemia (APL) cells, which suggests its participation in hematopoietic cell differentiation. In addition, ANXA8 is an anticoagulant protein that acts as an indirect inhibitor of the thromboplastin-specific complex, which is involved in the blood coagulation cascade.

配方:PBS

分子量:36.9 kDa

序列:Met 1-Pro327

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.