



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Zymogen Granule Membrane Protein 16/ZG16 Protein(C-6His)

产品货号: TP05435

第三版

别名:Zymogen Granule Membrane Protein 16; Zymogen Granule Protein 16; hZG16; Secretary Lectin ZG16; ZG16

描述:Zymogen Granule Membrane Protein 16 (ZG16) belongs to the jacalin lectin family. ZG16 is highly expressed in liver and is detected at lower levels in colon, ileum and jejunum. ZG16 may play a role in protein trafficking. In addition, ZG16 may act as a linker molecule between the submembranous matrix on the luminal side of zymogen granule membrane (ZGM) and aggregated secretory proteins during granule formation in the TGN.

配方:PBS

分子量:17.7 kDa

序列:Asn17-Cys167

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.