

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Interleukin-22 receptor subunit alpha-2/Interleukin-22 Binding Protein

产品货号: TP05440

第三版

**别名:** Interleukin-22 Receptor Subunit Alpha-2; IL-22 Receptor Subunit Alpha-2; IL-22R-Alpha-2; IL-22RA2; Cytokine Receptor Class-II Member 10; Cytokine Receptor Family 2 Member 10; CRF2-10; Cytokine Receptor Family Type 2 Soluble 1; CRF2-S1; Interleukin-22-Binding Protein; IL-22BP; IL22BP; ZcytoR16; IL22RA2

**描述:** Interleukin-22 Receptor Subunit  $\alpha$  -2 (IL22RA2) belongs to the type II cytokine receptor family. IL22RA2 is a secreted protein and contains three fibronectin type-III domains. IL22RA2 is widely expressed in many tissues. IL22RA2 functions as an IL22 antagonist and may be important in the regulation of inflammatory response. Three alternatively spliced transcript variants encoding distinct isoforms have been described.

**配方:** PBS

**分子量:** 25.7 kDa

**序列:** Thr 22-Pro231

**纯度:** > 95% by HPLC

**浓度:**

**内毒素:** <1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.