

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Protein Phosphatase 1 Regulatory Subunit 14A/PPP1R14A Protein(C-6His)

产品货号: TP05452

第三版

别名:Protein Phosphatase 1 Regulatory Subunit 14A; 17 kDa PKC-Potentiated Inhibitory Protein of PP1; Protein Kinase C-Potentiated Inhibitor Protein of 17 kDa; CPI-17; PPP1R14A; CPI17; PPP1INL

描述:Protein Phosphatase 1 Regulatory Subunit 14A (PPP1R14A) belongs to the PP1 inhibitor family. PPP1R14A is mapped to chromosome 19q13.13-q13.2. PPP1R14A binds directly to protein kinase C and casein kinase I. Meantime, PPP1R14A is a phosphorylation-dependent inhibitor of smooth muscle myosin phosphatase. PPP1R14A is the inhibitor of PPP1CA. When phosphorylated, PPP1R14A has over 1000-fold higher inhibitory activity, creating a molecular switch for regulating the phosphorylation status of PPP1CA substrates and smooth muscle contraction. In addition, inhibition of PPP1R14A also enhances contraction of smooth muscle in the absence of increment of intracellular Ca²⁺ concentration.

配方:PBS

分子量:17.8 kDa

序列:Met 1-Pro147

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.