



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Protein Phosphatase 1G/PP1MG Protein(C-6His)

产品货号: TP05457

第三版

别名:Protein Phosphatase 1G; Protein Phosphatase 1C; Protein Phosphatase 2C Isoform Gamma; PP2C-Gamma; Protein Phosphatase Magnesium-Dependent 1 Gamma; PPM1G; PPM1C

描述:Protein Phosphatase 1G (PP1MG) is a cytoplasmic protein that belongs to the PP2C family. PPM1G is widely expressed; it is abundant in testis, skeletal muscle, and heart. PPM1G is a negative regulator of cell stress response pathways. PPM1G is responsible for the dephosphorylation of Pre-mRNA splicing factors, an important factor for the formation of functional spliceosome. PPM1G also plays a role in regulating cell cycle progression.

配方:PBS

分子量:26.2 kDa

序列:Met317-Asp546

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.