

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Leukocyte Elastase Inhibitor/Serpin B1/SERPINB1 Protein(C-6His)

产品货号: TP05462

第三版

别名:Leukocyte elastase inhibitor;SERPINB1; Monocyte/neutrophil elastase inhibitor; M/NEI; Peptidase inhibitor 2; PI-2

**描述:**SERPINB1 is a member of the serpin family and Ov-serpin subfamily. As protease inhibitors, serpins have an array of functions including regulating blood coagulation, fibrinolysis, the complement pathway, angiogenesis, inflammation, tumor suppression, extracellular matrix remodeling, and cell motility. SERPINB1 regulates the activity of the neutrophil proteases elastase, cathepsin G, proteinase-3, chymase, chymotrypsin, and kallikrein-3. Reactive bond 1 of SerpinB1 is specific for reaction with chymotrypsin-like protease such as cathepsin G, chymotrypsin or chymase. Reactive bond 2 of SerpinB1 is specific for reaction with elastase-like protease such as neutrophyl elastase, proteinase-3, pancreatic elastase or PSA. In addition, SERPINB1 also functions as a potent intracellular inhibitor of granzyme H.

**配方:**PBS

**分子量:**43.8 kDa

**序列:**Met 1-Pro379

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.