

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Pulmonary Surfactant-Associated Protein D/PSP-D Protein(C-6His)

产品货号: TP05477

第三版

别名:Pulmonary Surfactant-Associated Protein D; PSP-D; SP-D; Collectin-7; Lung Surfactant Protein D; SFTPD; COLEC7; PSPD; SFTP4

**描述:**Surfactant Pulmonary-Associated Protein D (SP-D) is a 43 kDa member of the collectin family of innate immune modulators. Its principal components consist of a collagen-like region and a C-terminal carbohydrate recognition domain (CRD), a structure that places it in a subset of pattern recognition proteins termed defense collagens. SP-D is constitutively secreted by alveolar lining cells and epithelium associated with tubular structures and induced in cardiac smooth muscle and endothelial cells. It binds both secreted and transmembrane proteins that transduce its function. It binds human neutrophil defensins, modulating influenza anti-viral defense. It binds MD-2/LY96, a secreted protein that cooperates with Toll-like receptors (TLRs) in the response of macrophages to bacterial lipopolysaccharides (LPS) or cell wall components. It also binds macrophage CD14 and TLRs directly, blocking binding of LPS and down-regulating TNF- secretion. SP-D binding of both SIRP and the calreticulin/CD91 complex on macrophages allows for a graded response to environmental challenge.

**配方:**PBS

**分子量:**36.5 kDa

**序列:**Ala21-Phe375(Glu22Gly)?

**纯度:**> 95% by HPLC

**浓度:**

**内毒素:**<1.0 EU per 1 ug of protein (determined by LAL method)

**存储:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.