



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Brevican Core Protein/BCAN/BEHAB Protein(C-6His)

产品货号: TP05482

第三版

别名: Brevican Core Protein; Brain-Enriched Hyaluronan-Binding Protein; BEHAB; Chondroitin Sulfate Proteoglycan 7; BCAN; BEHAB; CSPG7

描述: Brevican Core Protein (BCAN) is a secreted protein that belongs to the aggrecan/versican proteoglycan family. BCAN contains one C-type lectin domain, one EGF-like domain, one Ig-like V-type domain, one Sushi (CCP/SCR) domain and two Link domains. BCAN may play a role in the terminally differentiating and the adult nervous system during postnatal development. BCAN could stabilize interactions between hyaluronan (HA) and brain proteoglycans.

配方: PBS

分子量: 97.8 kDa

序列: Asp23-Pro911

纯度: > 95% by HPLC

浓度:

内毒素: <1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.