

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Calnexin/CANX Protein(C-6His)

产品货号: TP05504

第三版

别名:Calnexin; IP90; Major Histocompatibility Complex Class I Antigen-Binding Protein p88; p90; CANX

描述:Calnexin/CANX is a single-pass type I membrane protein which belongs to the calreticulin family. It consists of a large N-terminal calcium-binding luminal domain, a single transmembrane helix and a short (90 residues), acidic cytoplasmic tail. The function of calnexin is to retain unfolded or unassembled N-linked glycoproteins in the endoplasmic reticulum. Calnexin is a calcium-binding protein that interacts briefly with newly synthesized glycoproteins in the endoplasmic reticulum. Calnexin may act in assisting protein assembly and/or in the retention within the ER of unassembled protein subunits. Calnexin seems to play a major role in the quality control apparatus of the ER by the retention of incorrectly folded proteins. Calnexin dwindles with aging and might contribute to a cytoprotection in an array of human age-related diseases.

配方:PBS

分子量:53.5 kDa

序列:His21-Pro481

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.