

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Follistatin-Related Protein 1/FSTL1 Protein(C-6His, E. coli)

产品货号: TP05516

第三版

别名:Follistatin-Related Protein 1; Follistatin-Like Protein 1; FSTL1; FRP

描述:Follistatin-Related Protein 1 (FSTL1) is a secreted protein that contains two EF-hand domains, one follistatin-like domain, one Kazal-like domain, and one VWFC domain. Its functional significance in physiological and pathological processes is not completely understood. However, FSTL1 is thought to modulate the action of some growth factors on cell proliferation and differentiation. FSTL1 maybe an autoantigen associated with rheumatoid arthritis.

配方:PBS

分子量:33.8 kDa

序列:Glu21-Ile308

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.