



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Serpin A10/ZPI Protein(C-6His)

产品货号: TP05525

第三版

别名:Protein Z-Dependent Protease Inhibitor; PZ-Dependent Protease Inhibitor; PZI; Serpin A10; SERPINA10; ZPI

描述:Serpin A10 is a secreted protein that belongs to the serpin family. It is predominantly expressed in the liver and secreted in the plasma. Its phosphorylation sites are present in the extracellular medium. It inhibits factors Xa and XIa of the coagulation cascade in the presence of protein Z, calcium, and phospholipid. Mutations in SERPINA10 are associated with venous thrombosis.

配方:PBS

分子量:49.5 kDa

序列:Leu22-Leu444

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.