

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Recombinant Human Serpin A6 /Transcortin Protein(C-6His)

产品货号: TP05527

第三版

别名:Corticosteroid-Binding Globulin; CBG; Serpin A6; Transcortin; SERPINA6; CBG

描述:Serpin Peptidase Inhibitor, Clade A, Member 6 (SerpinA6) belongs to the Serine Protease Inhibitors superfamily. SerpinA6 is synthesized in liver and has also been identified in a number of glucocorticoid responsive cells. SerpinA6 has an alpha-globulin protein with corticosteroid-binding properties. It is the major transport protein for glucocorticoids and progestins in the blood of most vertebrates. Defects in SERPINA6 are a cause of corticosteroid-binding globulin deficiency which is an extremely rare hereditary disorder characterized by reduced corticosteroid-binding capacity with normal or low plasma corticosteroid-binding globulin concentration, and normal or low basal cortisol levels associated with hypo/hypertension and muscle fatigue.

配方:PBS

分子量:43.7 kDa

序列:Met23-Val405

纯度:> 95% by HPLC

浓度:

内毒素:<1.0 EU per 1 ug of protein (determined by LAL method)

存储: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.